
HEAD OFFICE
Litron Lasers Ltd
8 Consul Road

Rugby Warwickshire CV21 1PB England

T +44 (0)1788 574444 F +44 (0)1788 574888 E sales@litron.co.uk

www.litronlasers.com





# LPY7000 High Energy Pulsed Nd:YAG Lasers Power, precision, performance.

#### **FEATURES**

- Rugged Industrial Build
- Up to 3.5J @ 1064nm
- Telescopic or Gaussian Resonators
- Optional Seeder Package
- All Harmonics to 5th Available
- Full RS232 Software Control
- Custom Systems up to 10J

#### **APPLICATIONS INCLUDE**

- Dye laser pumping
- OPO pumping
- Spectroscopy
- LIBS
- LIDAR
- PIV



The LPY7000 lasers offer extremely high
Q-switched outputs up to 3.5J and repetition
rates of up to 50Hz. Based around our proven
self-supporting invar frame their robust build
quality suits them to both industrial and
scientific applications. Energies up to 10J are
also available at 1064nm with extra
amplification.

The lasers are provided in an oscillator,


pre-amplifier, main-amplifier arrangement. The oscillator may be configured as a stable-telescopic resonator offering a low order multimode output

with a smooth spatial and temporal profile, or as an unstable Gaussian-coupled resonator offering a single transverse mode output with slightly higher peak powers.

Lamp change is performed in a matter of minutes with no need for any re-alignment at all. An optional IP54 sealed case ensures that the laser is protected against the ingress of dirt and moisture when used in industrial environments.



LPY7000 with additional 10J output stage doubled to 532nm outputting 5J at 1Hz. beam diameter 25mm.







# LPY7000 High Energy Pulsed Nd:YAG Lasers Power, precision, performance uncovered.

1 Rear Mirror


Intracavity Telescope
The intracavity telescope has a twofold use. Firstly it compensates for the thermal lensing of the laser rod. Secondly it reduces the intra-cavity beam diameter, thus effectively increasing the diffraction length in


the resonator. The result is very low divergence output beams whose beam profiles are spatially extremely homogeneous.

3 Electro-optic Q-Switch
A KD\*P Pockels cell is used within the Q-switch assembly.

Oscillator Pump Chamber
All pump chambers are machined from solid 316 stainless steel. The chambers are split such that the lamp housing can be removed easily during lamp changes, leaving the rod untouched. The chambers are fitted with close coupled ceramic reflectors for efficient and uniform pumping. A

proprietary filter plate protects the laser rod from UV emission from the lamp and also in the event of a lamp failure. The design of the chambers is such that a large turbulent water flow leads to very uniform cooling of the rod, essential for good pointing and overall stability.





An electronically verified, electronically actuated, intra-cavity safety shutter is standard on all of Litron's lasers.

6 Output Coupler

7 Expanding Telescope
An expanding telescope is used to expand and collimate the oscillator output prior to amplification.

8 Preamplifier Pump Chamber

9 Steering Mirrors

Expanding Telescope

Main Amplifier Pump Chambers

The main amplifier is configured in a birefringence-compensated twin-rod topology. This minimises the depolarisation of the laser beam and leads to more uniform and more efficient harmonic generation.

Invar Rail
The lasers are built on a rugged self supporting invar rail. This feature

sets them apart from all competitors as it is both more robust and more stable them conventional base-plate constructions. The modular nature of the rail allows for easy customisation of the lasers.

13 Harmonic Generation Unit

Output Stearing Mirrors

## LPY7000 High Energy Pulsed Nd:YAG Lasers

### **SPECIFICATIONS**

Model	LPY7864-10	LPY7864-20	LPY7864-30	LPY7864-50	LPY7875-10	LPY7875-20
Repetition Rate (Hz)	10	20	30	50	10	20
<b>Output Energy (mJ)</b> 1064nm 532nm 355nm 266nm	2750 1400 600 250	2250 1100 480 140	2000 1000 450 95	1400 700 250 80	3500 1750 700 275	2750 1350 600 170
Parameter Divergence (mrad) (1) Pointing Stability (urad) (2) Pulse length (ns) Jitter (ps) (3) Stability (+/-%) Beam diameter (mm) (4)  Services Voltage (VAC) Frequency (Hz) Water (I/min) @ 3-5bar	0.5 50 10-12 500 2 12.5 220-250 50-60 5	0.5 50 10-12 500 2 12.5 220-250 50-60 5	0.5 50 10-12 500 2 12.5	0.5 50 10-12 500 2 12.5 220-250 50-60 5	0.5 50 10-12 500 2 15 220-250 50-60 5	0.5 50 10-12 500 2 15 220-250 50-60 5

- Note
  (1) Full angle for 90% of the energy.
  (2) Full angle.
  (3) With respect to the external Q-switch trigger input.
  (4) Quoted as the main amplifier rod diameter.

### **MECHANICAL DATA**

